

CORONARY AND PERIPHERAL CATHETERS

| Name | Model Number | Minimum Vessel Diameter* | Compatibility | | Maximum OD Catheter | Working Length | Maximum Guide Wire |
|-------------|------------------------|--------------------------|---------------|-------|------------------------------|----------------|--------------------|
| | | | Sheath | Guide | | | |
| CP5 | CP100 (MN3006F) | 2.5–3.5 mm | 5F | 6F | 4.95 F (.065" or 1.65 mm) | 135 cm | .014" |
| CP55 | CP200 (MN3007F) | 3.5–5 mm | 6F | 7F | 5.33 F (.072" or 1.78 mm) | 135 cm | .014" |

PERIPHERAL CATHETERS

| | | | | | | | |
|---------------|---------------------|-----------|----|-----|---------------------------|--------|-------|
| PR7 | PR200 (5050) | 5–10.7 mm | 7F | n/a | 7 F (.093" or 2.36 mm) | 135 cm | .014" |
| PR7-65 | PR300 (5030) | 5–10.7 mm | 7F | n/a | 7 F (.093" or 2.36 mm) | 65 cm | .014" |

* Minimum vessel diameter references are recommendations