



It all comes out in the wash.

Rinspiration[®] may reduce the risk of peripheral occlusions pre- and post-intervention.

Peripheral Rinspirator

- Infusion holes at distal end direct the rinsing spray
- Central lumen aspirates debris
- Simultaneously rinses and aspirates thrombi from vessels
- Clears the vessel of mural thrombus, allowing better visibility of the lesion
- Delivers physician-specified solutions to the treatment site
- Provides better physician control with an easy-to-use, hand-held device
- Requires no capital equipment



U.S. Ordering Information for Coronary Rinspiration Catheter System

Model Number	Description	Catheter Length	Max. OD Catheter	Max. Guide Wire	Sheath Compatibility	Guide Compatibility
PR200	PR7 Peripheral Rinspirator 135cm	135cm	7F (0.093" or 2.36mm)	0.014"	7F	NA
PR300	PR7-65 Peripheral Rinspirator 65cm	65cm	7F (0.093" or 2.36mm)	0.014"	7F	NA
CP200	CP55 Coronary / Peripheral Rinspirator	135cm	5.33F (0.072" or 1.78 mm)	0.014"	6F	7F
CP100	CP5 Coronary / Peripheral Rinspirator	135cm	4.95F (0.065" or 1.65mm)	0.014"	5F	6F