



Introducing **SilverHawk[®] SXL**
In it for the long haul.

SilverHawk[®] SXL

Small Vessel Xtra Long Tip



SilverHawk SXL with maximum tissue capacity for small distal arteries (7.2cm tip)

- Increased tip length and micro vent holes maximize tissue capture
- Improved torque shaft allows for consistent torquing, tracking and device advancement
- Softer tip profile enhances maneuverability through distal vessels
- New Tissue Removal Device (TRD) provides efficient tip evacuation



This device incorporates MEC™ (Micro Efficient Compression) Technology - a breakthrough advancement that uses tiny vent holes to release fluid, and allows for greater tissue compression and storage.

SilverHawk[®] SXL

SilverHawk SXL

7.2cm Tip Length

P4033

Crossing Profile: 0.095" (2.4mm)

Guidewire Compatibility: 0.014"

Working Length: 136cm

Customer Service Phone: (877) FHT-HAWK

Customer Service Fax: (888) 233-2849

www.foxhollowtech.com

FOXHOLLOW

©2007 FoxHollow Technologies Inc.

LUS 422-A