


Clear the present danger.



Rinspiration[®] optimizes coronary intervention.

Coronary Rinspirator

- Infusion holes at distal end direct the rinsing spray
- Central lumen aspirates debris
- Simultaneously rinses and aspirates thrombi from vessels
- Clears the vessel of mural thrombus, allowing better visibility of the lesion
- Delivers physician-specified solutions to the treatment site
- Provides better physician control with an easy-to-use, hand-held device
- Requires no capital equipment



U.S. Ordering Information for Coronary Rinspiration Catheter System

Model Number	Description	Catheter Length	Max. OD Catheter	Max. Guide Wire	Sheath Compatibility	Guide Compatibility
CP100	CP5 Coronary / Peripheral Rinspirator	135cm	4.95F (0.065" or 1.65mm)	0.014"	5F	6F
CP200	CP55 Coronary / Peripheral Rinspirator	135cm	5.33F (0.072" or 1.78 mm)	0.014"	6F	7F